RESIST CTX-M

www.corisbio.com

IFU- 58R14/DE/05

Hersteller:

Coris BioConcept
CREALYS Science Park
Rue Guillaume Fouquet, 11
5032 GEMBLOUX
BELGIEN
Tel.: +32(0)81.719.917
Fax: +32(0)81.719.919
info@corisbio.com
Hergestellt in BELGIEN

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XII. PROBLEMES TECHNIQUES/RECLAMATIONS

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In vitro rapid diagnostic test for the detection of CTX-M-14 and CTX-M-15 β -lactamases in bacterial culture

FOR IN VITRO DIAGNOSTIC USE
FOR PROFESSIONAL USE ONLY

EN

References: K-15R14, 20 cassettes, buffer, 20 tubes and 20 transfer pipets

(EN) For Instructions For Use in your language : (FR) Pour obtenir les notices dans la langue de votre choix : (ES) Para las instrucciones de uso en su idioma : (PT) Para Instruções de Uso na sua língua : (IT) Per le Istruzioni di Uso nella sua lingua : (DE) Für Gebrauchsanleitungen in Ihrer Sprache : (NL) Voor Gebruiksaanwijzing in uw eigen taal : (PL) Instrukcja obsługi w języku użytkownika :	website: ifu.corisbio.com key-code: COR58R14
	(EU) +800 135 79 135 (non-EU) +31 20 794 7071 (CA) +1 855 805 8539 (AR, CO, UY, AU, NZ) +800 135 79 135

I. INTRODUCTION

β -lactamases are enzymes that inactivate β -lactams by cleaving the β -lactam ring, thereby inducing resistance to this type of antibiotic in bacteria producing these enzymes. Among the β -lactamases, those with an extended spectrum (ESBL) confer resistance to most beta-lactams, including 3rd generation cephalosporins, reducing the therapeutic window for the treatment of bacterial infections. The ESBLs of the CTX-M family belong to class A of β -lactamases, just like the enzymes TEM, SHV and the KPC type carbapenemases. Due to their evolution and particularly thanks to their association with mobile genetic elements, CTX-M enzymes represent an important part of ESBLs and are disseminated in many Enterobacteriaceae, in particular *Escherichia coli*, the genera *Klebsiella*, *Proteus*, *Citrobacter*, but also in other Enterobacteriaceae like *Morganella* spp.

There are over 250 variants of CTX-M described. They are classified into 3 minor and 5 major families among which group 1 (whose major representative is CTX-M-15) and group 9 (whose major representatives are CTX-M-14 and CTX-M-27) are the main responsible for the increase of ESBLs prevalence worldwide, whether in community or hospital settings. In contrast, the CTX-M-8-like cluster is a very small group of only 14 members.

Test for identification of CTX-M which represent the major ESBL part is a tool allowing a more rational use of last resort antibiotics such as carbapenems.

II. PRINCIPLE OF THE TEST

This test is ready to use and is based on a membrane technology with colloidal gold nanoparticles. Our kit is aimed at the detection of the group 1 (CTX-M-15-like), the group 8 (CTX-M-8-like) and the group 9 (CTX-M-14-like) CTX-M-type β -lactamases from a bacterial colony of Enterobacteriaceae growing on agar plate.

A nitrocellulose membrane is sensitized with (1) a monoclonal antibody directed against one epitope of the group 1 (CTX-M-15) β -lactamase and with (2) a second monoclonal antibody directed against a common epitope of group 8 (CTX-M-8-like) and group 9 (CTX-M-14-like) β -lactamases.

There are two different conjugates coupled to colloidal gold particles which are dried on a membrane: a conjugate directed against a second epitope of the CTX-M β -lactamase and a control conjugate to validate the test conditions.

When the provided buffer containing the resuspended bacteria comes into contact with the membrane, the solubilised conjugates migrate with the sample by passive diffusion and both conjugates and sample material come into contact with the first anti-CTX-M-15 antibody that it is adsorbed onto the nitrocellulose strip. If the sample contains a CTX-M-15 β -lactamase (or other CTX-M enzymes from group 1), the conjugate-CTX-M-15 complex will remain bound to the first (bottom) line anti-CTX-M-15 antibody adsorbed onto the nitrocellulose (line G1). The migration continues by passive diffusion and both conjugates and sample material come into contact with the second anti-CTX-M antibody line that is adsorbed onto the nitrocellulose strip. If the sample contains a

CTX-M-14 β -lactamase (or other CTX-M enzymes from group 9) or a CTX-M-8 β -lactamase, the conjugate-CTX-M complex will remain bound to the second line anti-CTX-M antibody adsorbed onto the nitrocellulose (line GX).

Finally, the solution continues to migrate to reach a third (upper) line control reagent that binds a control conjugate (line C), thereby producing a red line.

The result is visible within 15 minutes in the form of red lines on the strip.

III. REAGENTS AND MATERIALS

1. RESIST CTX-M (20)

20 sealed pouches containing one device and one desiccant. Each device contains one sensitized strip.

2. LY-E buffer vial (5 mL)

Saline solution containing Na₃N (<0,1%) and a detergent.

3. Instruction for use (1)

4. Disposable collection tubes (20)

5. Disposable transfer pipettes (20)

IV. SPECIAL PRECAUTIONS

- All operations linked to the use of the test must be performed in accordance with good laboratory practices.

- All reagents are for *in vitro* diagnostic use only.

- Pouch must be opened with care:

- Avoid touching nitrocellulose with your fingers.

- Wear gloves when handling samples.

- Never use reagents from another kit.

- Green lines indicate immunoreagents adsorption sites. Green colour disappears during the test.

- Reagents' quality cannot be guaranteed beyond their shelf-life dates or if reagents are not stored under required conditions as indicated in the insert.

V. WASTE DISPOSAL

- Dispose of gloves, swabs, test tubes and used devices in accordance with good laboratory practices.

- Each user is responsible for the management of any waste produced, and must ensure that it is disposed of in accordance with the applicable legislation.

VI. STORAGE

- An unopened pouch may be kept at between 4 and 30°C and used until the shelf-life date indicated on the packaging. Once the pouch is opened, perform the test immediately.

- Avoid freezing devices and buffer.

VII. SPECIMEN HANDLING AND COLLECTION

Specimens to be tested should be obtained and handled by standard microbiological methods.

Make sure that the specimens are not treated with solutions containing formaldehyde or its derivatives.

Culture media tested and validated with Coris BioConcept RESIT kits are listed on the website: <https://www.corisbio.com/faq/resist-ctx-m>

VIII. PROCEDURE

PREPARATIONS OF THE TEST:


Allow kit components, in unopened packaging, and specimens (in case the plate containing colonies to be tested was kept at 4°C) to reach room temperature (15-30°C) before performing a test.

Open the pouch and remove the device. Once opened, perform the test immediately. Indicate the patient's name or specimen number on the device (one device per sample).

SPECIMEN PREPARATION PROCEDURE:

Performance claims with regard to samples types other than bacterial colonies have not been established. We recommend the use of fresh bacterial colonies for optimal test performance.

1. Prepare one collection tube and add 4 drops of LY-E buffer in the tube.
2. Harvest bacteria by taking 3 colonies (equivalent to 3 loops of 1- μ L capacity) with a disposable bacteriological loop and dip the loop in the bottom of the collection tube containing the buffer.
3. Stir thoroughly before removing the loop
4. Close the tube and vortex the preparation to homogenize.
5. Use the transfer pipette provided in the kit and add 100 μ L of diluted sample into the sample well of the cassette (diluted sample must reach the black line indicated on the transfer pipette to accurately aspirate 100 μ L).
6. Allow to react for 15 minutes and read the result.

Positive results may be reported sooner the moment the test and control lines become visible.

Do not take the appearance of new lines into account after the reaction time is passed.

The result must be read on still wet strip.

Negative predictive value: 93.5 % (77.2 to 98.9 %)
Agreement: 98.6 % (146/148)

RESIST CTX-M	Groups X ^b	Positive	Negative	Total
Positive		26	0	26
Negative		0	122	122
Total		26	122	148

^b: Including CTX-M-8 and CTX-M-14

95 % Confidence Interval¹
Sensitivity: 100 % (83.4 to 100 %)
Specificity: 100 % (96.2 to 100 %)
Positive Predictive value: 100 % (84.0 to 100 %)
Negative predictive value: 100 % (96.2 to 100 %)
Agreement: 100 % (148/148)

C. Repeatability and reproducibility

To check intra-batch accuracy (repeatability), the same positive samples and a buffer solution were processed 15 times on kits of the same production batch in the same experimental conditions. All observed results were confirmed as expected. To check inter-batch accuracy (reproducibility), some samples (positive and buffer) were processed on kits from three different production batches. All results were confirmed as expected.

XI. LIMITS OF THE KIT

The test is qualitative and cannot predict the quantity of antigens present in the sample. Clinical presentation and other test results must be taken into consideration to establish diagnosis.

A positive test does not rule out the possibility that other antibiotic resistance mechanisms may be present.

XII. TECHNICAL PROBLEMS / COMPLAINTS

If you face a technical problem or if performances do not correspond with those indicated in this package insert:

1. Record the lot number of the kit concerned.
 2. If possible, keep the sample in the appropriate storage condition during the complaint management.
 3. Contact Coris BioConcept (client.care@corisbio.com) or your local distributor.
- Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the competent authority of the Member State in which the user and/or the patient is established.


XIII. BIBLIOGRAPHIC REFERENCES

- A. **Mendonça J., Guedes C., Silva C., Sá S., Oliveira M., Accioli G., Baylina P., Barata P., Pereira C., Fernandes R.** New CTX-M Group Conferring β -Lactam Resistance: A Compendium of Phylogenetic Insights from Biochemical, Molecular, and Structural Biology (Basel). 2022 Feb 7;11(2):256
- B. **Castanheira M, Simmer PJ, Bradford PA.** Extended-spectrum β -lactamases: an update on their characteristics, epidemiology and detection. JAC Antimicrob Resist. 2021 Jul 16;3(3):1-21
- C. **Peirano G, Pitout JDD.** Extended-Spectrum β -Lactamase-Producing Enterobacteriaceae: Update on Molecular Epidemiology and Treatment Options. Drugs. 2019 Sep;79(14):1529-1541
- D. **Rossolini GM, D'Andrea MM, Mugnaioli C.** The spread of CTX-M-type extended-spectrum beta-lactamases. Clin Microbiol Infect. 2008 Jan;14 Suppl 1:33-41
- E. **Ghafourian S, Sadeghifard N, Soheili S, Sekawi Z.** Extended Spectrum Beta-lactamases: Definition, Classification and Epidemiology. Curr Issues Mol Biol. 2015;17:11-21.
- F. **Paterson DL, Bonomo RA.** Extended-spectrum beta-lactamases: a clinical update. Clin Microbiol Rev. 2005 Oct;18(4):657-86
- G. **Bevan E.R., Jones A.M. and Hawkey P.M.** Global epidemiology of CTX-M β -lactamases: temporal and geographical shifts in genotype. J Antimicrob Chemother. 2017 Aug 1;72 (8):2145-2155.
- H. **Zhao WH, Hu ZQ.** Epidemiology and genetics of CTX-M extended-spectrum β -lactamases in Gram-negative bacteria. Crit Rev Microbiol. 2013 Feb;39(1):79-101.
- I. **D'Andrea MM, Arena F, Pallecchi L, Rossolini GM.** CTX-M-type β -lactamases: a successful story of antibiotic resistance. Int J Med Microbiol. 2013 Aug;303(6-7):305-17
- J. **Newcombe, Robert G.** "Two-Sided Confidence Intervals for the Single Proportion: Comparison of Seven Methods," *Statistics in Medicine*, 17, 857-872 (1998)

Last update : 27 MAY 2024

REF	Catalogue number		Manufacturer
IVD	In vitro diagnostic medical device		Temperature limits
	Contains sufficient for <n> tests	LOT	Batch code
	Consult instructions for use		Do not reuse
	Keep dry		Use by
DIL SPE	Diluent specimen	CONT NaN₃	Contains Sodium azide
UDI	Unique device identifier		

IX. INTERPRETING RESULTS

The results are to be interpreted as follows:

Negative test result: a reddish-purple line appears across the central reading window at the Control line (C) position. No other band is present.

Positive test result: in addition to a reddish-purple band at the Control line (C), a visible reddish-purple band appears at one of the Test line position ("G1" or "GX"). Intensity of the test line may vary according to the quantity of antigens as well as of the variant type present in the sample. Any reddish-purple Test line ("G1" or "GX"), even weak, should be considered as a positive result.

If only the positive test line is line GX, the sample contains CTX-M-14 β -lactamase or another CTX-M enzymes from group 9 or from group 8. If only the positive test line is line G1, the sample contains CTX-M-15 β -lactamase or another CTX-M enzymes from group 1.

Invalid test result: The absence of a Control line indicates a failure in the test procedure. Repeat invalid tests with a new test device.

Note: during the drying process, a very faint shadow may appear at the Test line positions. It should not be regarded as a positive result.

X. PERFORMANCE

A. Detection Limit

The detection limit determined with purified recombinant proteins of CTX-M-14 and CTX-M-15 have been evaluated at 3.14 ng/ml and 3.7 ng/ml respectively.

B. Validation on collection of reference strains

RESIST CTX-M test was evaluated on a collection of 148 clinical isolates with fully characterized beta-lactam antibiotic resistance mechanisms by a molecular approach (in house multiplex PCR and sequencing).

RESIST CTX-M	Group 1 ^a	Positive	Negative	Total
Positive		117	0	117
Negative		2	29	31
Total		119	29	148

^a:including CTX-M-1, 3, 15, 55, 135 and CTX-M G1 not sequenced

95 % Confidence Interval¹
Sensitivity: 98.3 % (93.6 to 99.7 %)
Specificity: 100 % (85.4 to 100 %)
Positive Predictive value: 100 % (96.0 to 100 %)